

Fixed height sampling chair

*512 range* **512 02T** 







#### **TECHNICAL DETAILS**

- 3 sections sampling chair with fixed height 67cm
- Head rest adjustment in every positions with gas jack
- Legrest adjustment via a self locking mechanism system, the legrest maybe positioned vertically
- Trendelenburg position
- Flat upholstery with good density 35kg/cm², M1 fire retardant class, thickness 8cm
- Equipped with 2 blood test splints, adjustable in every positions with serrated blocks system
- Epoxy covering, light grey colour (RAL 7035) in standard

### **DIMENSIONS**

- Width 55cm
- Height 67cm
- Load capacity: 130kg


## **OTHER ELEMENTS**

- All our chairs can be put in **lying position**
- Low-height sampling chairs (50cm) are intended for blood sampling, higher chairs are dedicated to gynecological sampling.

To give a better confort to patients and practitioners

### You may also like:

Ref. **512 02TB** : Same as 512 02T but height 50cm





# WHAT YOU NEED TO KNOW

- Delivery time: 15 days upon receipt of order
- Warranty: 2 years (1 year for the upholstery)

V01-10/21 Page 2 sur 7

#### **ACCESSORIES**



**ACC 537** 

Pair of clamps diam. 14mm for blood test splints and stirrups for all ranges except LUVIA



ACC 535

Pair of double clamps diam. 14 and diam. 16mmto fix both armrest and infusion stand





**ACC 964** 

Pair of legrests on knee system in black polyurethane, diam. 14 mm, without clamp



ACC 105

Pair of stainless steel stirrups diam. 14mm without clamp



Plastic protection



**SLFEMIREPOSE** 

Cover for legrests ACC964 in upholstery, the unit



ACC 514-GOUT

for blood test splint, the unit



ACC 514-ASSISE

Plastic protection for seatrest



**ACC 514-JAMBIÈRE** 

Plastic protection for legrest



ACC 517-512

Plastic protection for backrest



ACC 54005

Half moon head rest width 60cm fire retardant M1 class



ACC 910-PLAT02

Stainless steel tray 340x340 with clamp



**ACC 910-PLAT04** 

Plastic tray 350x270 with clamp

Pair of handles with clamps



**ACC 113U** 

Pair of handles without clamp with U shape



ACC 113U-105

to fix both handles and stirrups



ACC 125

Paper roll holder for all models



ACC 32012

Stainless steel infusion stand diam. 16mm, 2 plastic hooks, without clamp



ACC 32014

Stainless steel infusion stand diam. 16mm, 2 stainless steel hooks without clamp



ACC 9505

Set of 4 wheels with brake diam. 50mm



ACC 993-FIXE

4 retractable wheels diam. 50mm



**ACC TETERGO** 

Ergonomic headrest



ACC 110

Blood test splint in colour without clamp, the unit

V01-10/21 Page 3 sur 7

## UPHOLSTERY'S COLOURS

Fire retardant upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range Be careful! Urban colours **with extra value** 

## Upholstery: Classic range



## Upholstery: Urban range



V01-10/21 Page 4 sur 7




The entire 512 range is offered in light grey (RAL 7035) as standard Possibility to change the color without extra charge

#### Base colour



V01-10/21 Page 5 sur 7



#### **DAILY MAINTENANCE**

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

#### WARRANTY

Mechanical/electrical part: 2 years

Upholstery: 1 year

V01-10/21 Page 6 sur 7

## **CONTACT**



Parc d'entreprises Visionis 278 rue de l'industrie 01090 GUÉREINS- FRANCE

#### **COMMERCIAL DEPARTMENT**

Tél: +33(0)4 74 66 69 49 Fax: +33(0)4 74 66 69 50 export@carinamedical.com

V01-10/21 Page 7 sur 7