



Folding screen
116 range
116 01-PB

CARINA
Medical




CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- Folding screens with 3 panels
- Panels mounted on hinges
- Screen in white plastic, fire retardant M1 class
- Structure in square steel tubes with epoxy covering, silver or white color (RAL 9003)
- Equipped with 6 wheels diam. 50mm

DIMENSIONS

- Panels width 700x600x700mm
- Unfold length 2m
- Height 1,8m



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 3 years (1 year for the upholstery)

You may also like:

Réf. 114 01-PB: Same as 116 01-PB but with 2 panels



**BASE
COLOURS**



RAL 9003



RAL 1013



RAL 1015



RAL 6019



RAL 6033



RAL 6000




RAL 5017



RAL 5002



RAL 5000



RAL 7035



RAL 7016



RAL 9005




RAL 1018



RAL 3020



RAL 3003



SILVER



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

3 years

Upholstery: 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com