



Folding screens
Range 119
119 03-PB

CARINA
Medical





CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- Folding screens 3 panels
- Panels mounted on hinges
- White PVC, fire retardant, M1 class
- Structure in square steel tubes epoxy covering, silver or white color (RAL 9003)
- Equipped with 6 wheels diam. 50mm

DIMENSIONS

- Panels width 700x600x700mm
- Unfold length 2m
- Height 1,8m



WHAT YOU NEED TO KNOW

- ✓ Warranty : 3 years
(1 year for the upholstery)

You may also like :

Réf. **119 02-PB** : Same as 119 03-PB but 2 panels

Réf. **119 04-PB** : Same as 119 03-PB but 4 panels



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

3 years



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL DEPARTMENT

Tél : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com