



Stool
120 range
12001CHR

CARINA
Medical



CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- Height adjustable stool with chromed gas column
- Manual control under the seat
- Fire retardant upholstery M1 class, thickness 6cm
- Equipped with 5 wheels, diameter 50mm
- 5 branches chromed base diameter 550mm

DIMENSIONS

- Round seat, diameter 35cm
- Adjustable height from 48 to 68cm
- Load capacity : 130kg

You might also like :

Ref. **12011CHR** : Same as 12001CHR but on skids, variable height from 44 to 64cm



WHAT YOU NEED TO KNOW

- ✓ Delivery time : 15 days upon receipt of order
- ✓ Warranty : 2 years (1 year for upholstery)

ACCESSORIES



ACC 10
Backrest for PRISCA,
range 120, 30 and 40



ACC 14
Foot control



ACC 15
Footrest



ACC 16
Chromed base,
diameter 38cm
(Reducing height of 1cm)



ACCROULETTE-02
5 single-wheels,
diameter 50mm



ACCROULETTE-03
5 double-wheels for hard
floors, diameter 50mm



ACCROULETTE-04
5 antistatic wheels, diameter
50mm



ACCROULETTE-05
5 auto-braking wheels,
in load or out of load

HH/HB
High or low height

With ACC 10
option




UPHOLSTERY'S COLOURS

Fire retardant upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range

Upholstery : Classic range




Arno




Garonne




Potomac




Amour




Mississippi



St Laurent




Rhône



Missouri



Seine



Amazone

Upholstery : Urban Range



Abysses




Cendre



Origan




Géranium




Réglisse



Prune



Perle




Ficelle




Terre




Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate



MAINTENANCE GUIDE

DAILY MAINTENANCE

- Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvaniol pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material
- Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

- Mechanical / electrical part : 2 years
- Upholstery : 1 year



CONTACT

CARINA
Medical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tél : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com