



Stool
65 range
6506

CARINA
Medical



CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- Height adjustable stool with chromed gas column
- Manual control under the seat
- Saddle seat
- Fire retardant upholstery M1 class, thickness 6cm
- Equipped with 5 wheels, diameter 50mm
- 5 branches aluminium base diameter 500mm

DIMENSIONS

- Adjustable height from 50 to 70cm
- Load capacity : 180kg

You might also like :

Ref. **6507** : Same as 6506 but on skids, adjustable height from 46 to 66cm



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 2 years (1 year for upholstery)

ACCESSORIES



ACC 13B
Sanglar backrest
for 65 range

ACC 17
Front/rear seat tilt for 65 range



ACC 14
Foot control



ACC 15
Footrest



ACC 16
Chromed base,
diameter 38cm
(Reducing height of 1cm)



ACCCROULETTE-02
5 single-wheels,
diameter 50mm



ACCCROULETTE-03
5 double-wheels for hard
floors, diameter 50mm



ACCCROULETTE-04
5 antistatic wheels, diameter
50mm



ACCCROULETTE-05
5 auto-braking wheels,
in load or out of load

HH/HB
High or low height

with ACC 13B
option



UPHOLSTERY'S COLOURS

Fire retardant upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range

Upholstery : Classic range



Upholstery : Urban Range





MAINTENANCE GUIDE

DAILY MAINTENANCE

- Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvavios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material
- Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

- Mechanical / electrical part : 2 years
- Upholstery : 1 year



CONTACT

CARINA *Medical*

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tél : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com