



*Picture of the model with all around controller
and removable wheels (Ref. 99926)*

Physiotherapy table

999 range

999 22

CARINA
Medical



CE


ISO 9001
ISO 13485

TECHNICAL DETAILS

- 2 sections physiotherapy table, with electric height adjustment
- Pedal control
- Adjustment of the backrest in every position with a gas jack
- Stretched upholstery with good density 40g/cm², M2 fire retardant class, thickness 7 cm
- Height adjustment with an electric jack of 6000N capacity
- Face hole to give a better confort to patients, with plug cushion
- Equiped with exam sheet roller on the foot side
- Epoxy covered base, silver color in standard
- With 4 fixed feets

DIMENSIONS

- Width 600mm
- Adjustable height from 48 to 95cm
- Load capacity: 150kg



You may also like :

Réf. 99921 : Same as 99922 but with hand remote control

Réf. 99923 : Same as 99922 but with all around controller

Réf. 99924 : Same as 99922 but with hand remote control and 4 removable wheels

Réf. 99925 : Same as 99922 but with 4 removable wheels

Réf. 99926 : Same as 99922 but with all around controller and 4 removable wheels




WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 2 years (1 year for the upholstery)

ACCESSORIES



ACC 538
Pair of clamps diam.14mm
width 600mm



ACC 105
Pair of stainless steel stirrup
diam. 14mm,without clamp



ACC964
Pair of legrests on knee system
in black polyurethane,
diam.14mm, without clamp



SLFEMIREPOSE
Cover for legrests ACC964
in upholstery, the unit



ACC 76C
Unit clamp
to fix infusion stand diam.16mm




ACC32012
Stainless steel infusion stand
diam. 16mm
2 plastic hooks, without clamp



ACC 110
Blood test splint in colour
without clamp, the unit



ACC 110POLY
Polyurethane blood test splint
without clamp, the unit



ACC 54001
Half moon headrest width 600mm
fire retardant M1 class



ACCTETERGO
Ergonomic head rest



ACC 514-UNIVERSE
Plastic protection cover

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M2 class

For the upholstery colours, you can choose between the Ecomax, Classic and Urban range
Be careful ! Classic and Urban colours **with extra value**

Upholstery: ECOMAX* range



Galapagos



Bahamas



Blanc



Tahiti



Malouines


*Standard colors for ECOMAX, SIMPLEX and 999 ranges,
Classic and Urban colors available **on demand**

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range

Upholstery: Classic range



Arno




Garonne



Potomac




Amour




Mississippi




St Laurent



Rhône



Missouri



Seine




Amazone

Upholstery: Urban range



Abysses




Cendre




Origan




Géranium




Réglice




Prune




Perle




Ficelle




Terre




Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvaniol pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part: 2 years

Upholstery: 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com