



Fixed height sampling chair

522 range

522+990

CARINA
Medical




CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- 3 sections sampling chair with fixed height 50cm
- Backrest and legrest adjustment with a synchronised gas jack
- The legrest can be put in a vertical position
- Trendelenburg position
- Ergonomic upholstery with good density 35kg/cm², M1 fire retardant class, thickness 7cm. Specially developed for optimal patient comfort
- Equipped with 2 polyurethane test splints adjustable in every position with serrated blocks system and a push bar
- Chair with infusion stand holder and a pair of clamps to fix stirrups (option)
- Epoxy covering, light grey colour (RAL 7035) in standard
- On 4 wheels with brake diam. 125mm

DIMENSIONS

- Width 55cm
- Height 50cm
- Overall dimensions :
 - Width 80cm
 - Length when chair upright 88cm
- Load capacity : 180kg



OTHER ELEMENTS

All our chairs can be put in **lying position**

You may also like :

Ref. **522** : Same as 522+990 but without wheels



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty : 2 years (1 year for the upholstery)

ZOOM 522+990



ACCESSORIES



ACC 535
Pair of double clamps
diam. 14 and diam. 16mm to fix
both armrest and infusion stand

ACC 516

Option synchronized
backrest and legrest




ACC 964
Pair of legrests on knee system
in black polyurethane,
diam. 14 mm, without clamp



ACC 514-GOUT
Plastic protection
for blood test splint, the unit



ACC 514-522-JAMB
Plastic protection for legrest



ACC 54005
Half moon head rest width
60cm fire retardant M1 class

ACC 910-PLAT04

Plastic tray 350x270 with clamp



ACC 113U-105
Pair of handles with clamps
to fix both handles and
stirrups



ACC 32012
Stainless steel infusion stand
diam. 16mm, 2 plastic hooks,
without clamp



ACC 110
Blood test splint in colour
without clamp, the unit



ACC 125
Paper roll holder for all models



ACC 105
Pair of stainless steel
stirrups diam. 14mm
without clamp



SLFEMIREPOSE
Cover for legrests ACC964
in upholstery, the unit



ACC 514-ASSISE
Plastic protection for seatrest



ACC 517-512
Plastic protection for backrest



ACC 910-PLAT02
Stainless steel tray 340x340
with clamp



ACC 113U
Pair of handles without clamp
with U shape



ACC TETERGO
Ergonomic headrest



ACC 32014
Stainless steel infusion stand
diam. 16mm, 2 stainless steel
hooks without clamp



ACC 110POLY
Polyurethane blood test splint
without clamp, the unit

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range
Be careful ! Urban colours **with extra value**

Upholstery : Classic range




Arno




Garonne




Potomac



Amour



Mississippi



St Laurent




Rhône



Missouri



Seine



Amazone

Upholstery : Urban range



Abysses




Cendre



Origan




Géranium




Réglisse




Prune



Perle




Ficelle




Terre




Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate




BASE COLOURS

The entire 522 range is offered in light grey (RAL 7035) as standard
Possibility to change the color without extra charge


Base colour



RAL 9003



RAL 1013



RAL 1015



RAL 6019



RAL 6033



RAL 6000




RAL 5017



RAL 5002




RAL 5000



RAL 7035



RAL 7016



RAL 9005




RAL 1018



RAL 3020



RAL 3003



Shiny

SILVER



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part : 2 years

Upholstery : 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL DEPARTMENT

Tél : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com