



Fixed height examination couch

ECOMAX range

ECOMAX-700

CARINA
Medical




CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- 2 sections examination couch, width 70cm
- Adjustable backrest with 2 indexable racks
- Flat upholstery with good density 24kg/m², M1 fire retardant class, thickness 6cm
- Equipped with paper roll holder : possibility to fix it on foot or head side
- Epoxy covering, silver colour


DIMENSIONS

- Width 70cm
- Height 80cm
- Load capacity : 150kg



You may also like :

Ref. **ECOMAX** : Same as ECOMAX-700 but width 60cm



WHAT YOU NEED TO KNOW

- ✓ Delivered disassembled (assembly with extra value)
- ✓ Warranty : 5 years (1 year for the upholstery)

ACCESSORIES



ACC 538D16
Pair of clamps diam. 16mm for
ECOMAX range width 60cm



ACC 538-700
Pair of clamps diam. 14mm
for ECOMAX range width 70cm



ACC 76C
Unit clamp to fix infusion stand
diam. 16mm




ACC 32012
Stainless steel infusion stand
diam. 16mm, 2 plastic hooks,
without clamp



ACC 110
Blood test splint in colour
without clamp, the unit



ACC 110POLY
Polyurethane blood test splint
without clamp, the unit



ACC 105
Pair of stainless steel
stirrups diam. 14mm,
without clamp



ACC 964
Pair of legrests on knee system
in black polyurethane
diam. 14mm, without clamp




SLFEMIREPOSE
Cover for legrests (ACC964)
in upholstery, the unit



ACC 514-UNIVERSE
Plastic protection cover
for all couches



ACC TETERGO
Ergonomic head rest



ACC 54002
Half moon head rest width
70cm fire retardant M1 class



ACC 832
Face hole with plug cushion



ACC URGENCE03-70
Pair of retractable security side
rails for ECOMAX couch
width 70cm



ACC 990
Set of 4 wheels with brake
diam. 125mm

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M2 class

For the upholstery colours, you can choose between the Ecomax, Classic and Urban range
Be careful ! Classic and Urban colours **with extra value**


Upholstery : ECOMAX range*



Galapagos



Bahamas



Blanc



Tahiti



Malouines

*Standard colors for ECOMAX and SIMPLEX range,
Classic and Urban colors available **on demand**



UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range
Be careful ! Classic and Urban colours **with extra value**

Upholstery : Classic range




Arno



Garonne




Potomac



Amour




Mississippi



St Laurent




Rhône



Missouri



Seine



Amazone

Upholstery : Urban range



Abysses



Cendre



Origan




Géranium




Régilisse




Prune




Perle




Ficelle




Terre




Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part : 5 years

Upholstery : 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL DEPARTMENT

Tél : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com