

Hydraulic sampling chair

LUVIA range LUVIA14T








TECHNICAL DETAILS

- 3 sections hydraulic sampling chair, 360° rotating (possibility to block the rotation)
- Backrest and legrest adjustment via independent gas jacks
- The legrest can be put in a vertical position
- Trendelenburg position
- Ergonomic upholstery with good density 35kg/cm², M1 fire retardant class, thickness
 7cm. Specially developed for optimal patient comfort
- Equipped with 2 blood test splints, adjustable in every positions with serrated blocks system
- Plastic protection on seat and legrest as standard
- Chromed base

DIMENSIONS

- Width 55cm
- Adjustable height from 54cm to 73cm
- Load capacity: 130kg in dynamic





OTHER ELEMENTS

All our chairs can be put in **lying position**

You may also like:

Ref. **LUVIA14TSJ**: Same as LUVIA14T but

without legrest



WHAT YOU NEED TO KNOW

- Delivery time: 15 days upon receipt of order
- Warranty: 2 years (1 year for the upholstery)

V01-10/21 Page 2 sur 6

ACCESSORIES



ACC 537LUVIA

Pair of clamps diam. 14mm for blood test splints and stirrups for LUVIA range



ACC 105

Pair of stainless steel stirrups diam. 14mm, without clamp



SLFEMIREPOSE

Cover for legrests (ACC964) in upholstery, the unit



ACC514-LUVI-ASS

Plastic protection for seatrest



ACC514-LUVI-DOS

Plastic protection for backrest



ACC 910-PLAT02

Stainless steel tray 340x340 with clamp



ACC 113U

Pair of handles without clamp with U shape



ACC 32012

Stainless steel infusion stand diam. 16mm, 2 plastic hooks, without clamp



ACC TETERGO-LUV

Ergonomic head rest for LUVIA range



ACC 110POLY

Polyurethane blood test splint without clamp, the unit



ACC 535

Pair of double clamps diam. 14mm and diam.16mm to fix both armrest and infusion stand



ACC 964

Pair of legrests on knee system in black polyurethane diam. 14mm, without clamp



ACC 514-GOUT

Plastic protection for blood test splint, the unit




ACC514-LUVI-JAM

Plastic protection for legrest



ACC 54005

Half moon head rest width 55cm fire retardant M1 class



Plastic tray 350x270 with clamp



ACC 113U-105

Pair of handles with clamps to fix both handles and stirrups




ACC 32014

Stainless steel infusion stand diam. 16mm, 2 stainless steel hooks without clamp



ACC 110

Blood test splint in colour without clamp, the unit



ACC 125

Paper holder for all models

V01-10/21 Page 3 sur 6

UPHOLSTERY'S COLOURS

Fire retardant upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range Be careful! Urban colours **with extra value**

Upholstery: Classic range



Upholstery: Urban range



V01-10/21 Page 4 sur 6



DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/electrical part: 2 years

Upholstery:1 year

V01-10/21 Page 5 sur 6





Parc d'entreprises Visionis 278 rue de l'industrie 01090 GUÉREINS- FRANCE

COMMERCIAL DEPARTMENT

Tél: +33(0)4 74 66 69 49 Fax: +33(0)4 74 66 69 50 export@carinamedical.com

V01-10/21 Page 6 sur 6