



Fixed height sampling chair

PRELIA range

PRELIA07

CARINA
Medical




CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- 2 sections sampling chair with fixed height 50cm
 - Adjustment of the flexibility of the backrest
 - Ergonomic upholstery, M1 fire retardant class, thickness 6cm.
- Very good lumbar support**
- Equipped with a blood test splint in black polyurethane
 - 5 aluminium branches base diam. 65 cm on skids

DIMENSIONS

- Width 46cm
- Height 50cm
- Load capacity : 150kg



ACCESSORIES



ACC 110POLY
Polyurethane blood test splint
without clamp, the unit



ACC 110
Blood test splint in colour
without clamp, the unit



ACC 537-PRELI-UN
Pair of clamps diam. 14mm for
blood test splints and stirrups
for PRELIA range

You may also like :

Ref. **PRELIA08** : Same as PRELIA07 but
equipped with a blood test splint, M1 fire
retardant class



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days
upon receipt of order
- ✓ Warranty : 2 years (1 year
for the upholstery)

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range
Be careful ! Urban colours without extra value

Upholstery : Classic range



Arno




Garonne




Potomac




Amour




Mississippi



St Laurent




Rhône



Missouri



Seine



Amazone



Discover the Urban range on the
following page

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

Discover the upholstery colors of the Urban range



Upholstery : Urban range





MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part : 2 years

Upholstery : 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL DEPARTMENT

Tél : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com