



*Electric examination couch
SIMPLEX range*

SIMPLEX-3.0


CARINA
Medical



CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- 2 sections examination couch, width 60cm
- Control: Pedal (hand remote on request)
- Backrest adjustment in every positions with gas jack
- Flat upholstery with good density 24kg/cm², M2 fire retardant class, thickness 6cm
- Equipped with paper roll holder for exam sheet, on the foot side
- Epoxy covering, silver colour



DIMENSIONS

- Width 60cm
- Adjustable height from 48 to 88cm
- Load capacity: 200kg in dynamic

Couch available with wheels
option ACC990-SIMPLEX
Ref: SIMPLEX-.3.0+990



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 2 years (1 year for the upholstery)

You may also like :

Ref. **SIMPLEX-3.0-700**: Same as SIMPLEX-3.0 but width 70cm

ACCESSORIES



ACC 990-SIMPLEX

Set of 4 wheels with brake
diam 75mm




ACC 76C

Unit clamp to fix infusion stand
diam 16mm



ACC 110

Blood test splint in colour
without clamp, the unit



ACC 105

Pair of stainless steel stirrups
diam 14mm, without clamp



SLFEMIREPOSE

Cover for legrests (ACC964) in
upholstery, the unit



ACC 54001

Half moon head rest width 60cm
fire retardant M1 class



ACC URGENCE03-60

Pair of retractable security
side rails for SIMPLEX couch
width 60cm



ACC 832

Face hole with plug cushion



ACC 538

Pair of clamps diam 14mm for
SIMPLEX range width 60cm



ACC 32012

Stainless steel infusion stand
diam 16mm, 2 plastic hooks
without clamp



ACC 110POLY

Polyurethane blood test splint
without clamp, the unit



ACC 964

Pair of legrests on knee system
in black polyurethane diam
14mm, without clamp



ACC 514-UNIVERSE

Plastic protection cover
for all couches



ACC TETERGO

Ergonomic headrest




ACC BARRILITCOM

Pair of retractable
security side rails



BASE COLOURS



Fire retardant
upholstery, M2 class

*For the upholstery colours, you can choose between the Ecomax, Classic and Urban range
Be careful ! Classic and Urban colours with extra value*

Upholstery: ECOMAX range*



Galapagos



Bahamas



Blanc



Tahiti



Malouines


**Standard colors for ECOMAX and SIMPLEX range,
Classic and Urban colors available on demand*

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

*For the upholstery colours, you can choose between the Classic and Urban range
Be careful ! Urban colours with extra value*

Upholstery: Classic range



Arno



Garonne



Potomac



Amour



Mississippi




St Laurent



Rhône



Missouri



Seine



Amazon

Upholstery: Urban range



Abyss



Cendre



Origan




Géranium




Régisse




Prune




Perle



Ficelle



Terre




Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part: 2 years

Upholstery: 1 year



CONTACT



Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com