



Picture of the model with 3 panels

Folding screen

118 range

118 04

CARINA
Medical



CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- Folding screens with 4 panels
- Panels assembled with zip
- Fire retardant upholstery, M1 class
- Structure in round steel tubes with epoxy covering, light grey color (RAL 7035) in standard
- Equipped with 13 wheels diam. 50mm

DIMENSIONS

- Panels width 700x600x600x700mm
- Unfold length 2,60m
- Height 1,8m



**Panels can be unassembled
for an easy cleaning**

WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 3 years (1 year for the upholstery)

You may also like:

Réf. 118 03: Same as 118 04 but with 3 panels


Réf. 118 02: Same as 118 04 but with 2 panels

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

*For the upholstery colours, you can choose between the Classic and Urban range
Be careful ! Urban colours with extra value*

Upholstery: Classic range



Arno




Garonne




Potomac




Amour



Mississippi



St Laurent




Rhône



Missouri



Seine




Amazone

Upholstery: Urban range



Abysses



Cendre



Origan




Géranium




Réglisse




Prune




Perle




Ficelle




Terre



Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate



**BASE
COLOURS**




RAL 9003



RAL 1013



RAL 1015



RAL 6019



RAL 6033



RAL 6000




RAL 5017



RAL 5002



RAL 5000



RAL 7035



RAL 7016



RAL 9005




RAL 1018



RAL 3020



RAL 3003



SILVER



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvaniol pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

3 years

Upholstery: 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com