



*Infrared lamp 250W*

**172 range**

**172 01**

**CARINA**  
*Médical*



**CE**  
ISO 9001  
ISO 13485

## TECHNICAL DETAILS

- Halogen lamp 250 Watts, white colour
- Flexible arm, length 860mm (including 340mm rigid)
- Wavelength: 1200nm
- Infrared bulb 250W ( supplied in standard) with an E27 base and integrated reflector: transmits 90% of the energy in infrared
- Bulb duration: 5000h minimum
- Ergonomic handle to easy positioning of the lamp
- Head diam. 260mm mounted on orientation joint
- Power cable length : 4m
- Epoxy covered base, white colour (RAL 9003)
- 5 rotating wheels
- Weight: 3kg
- Classe I medical device

### OPTION

ISC (no contact switch- Timer from 0 to 60 minutes)

**Ideal lamp  
physiotherapy,  
osteopathy, chiropractic**



### WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 5 years



CONTACT



Parc d'entreprises Visionis  
278 rue de l'industrie  
01090 GUÉREINS- FRANCE

**COMMERCIAL SERVICE**

Tel: +33(0)4 74 66 69 49  
Fax: +33(0)4 74 66 69 50  
[export@carinamedical.com](mailto:export@carinamedical.com)