



*Halogen lamp 20W*

**173 range**

**173 02**

**CARINA**  
*Medical*



**CE**  
ISO 9001  
ISO 13485

## TECHNICAL DETAILS

- Halogen lamp 20 Watts, black colour
- Flexible arm, length 650mm
- Light intensity: 300 Lumens-3000° Kelvin (neutral white)
- CRI: 100/100
- Luminous power: 23000 Lux/ 50cm, illuminated diameter: 100mm
- Consumption: 20W
- Dichroic source 12V
- Bulb longevity: 5000h minimum
- Head diam. 70mm mounted on orientation joint
- Power cable length : 3m
- Epoxy covered base, white or black colour (RAL 9003 or RAL 9005)
- 5 rotating wheels
- Classe I medical device

### OPTION

ISC (No contact switch)

**Ideal lamp  
for general medicine  
& nursing**

### You may also like :

Réf. **173 01**: Same as 173 01 but white colour

Réf. **173 03**: Same as 173 01 but wall-mounted

Réf. **173 04**: Same as 173 01 but white colour and wall-mounted



### WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 5 years



CONTACT

CARINA  
*Médical*

Parc d'entreprises Visionis  
278 rue de l'industrie  
01090 GUÉREINS- FRANCE

**COMMERCIAL SERVICE**

Tel: +33(0)4 74 66 69 49  
Fax: +33(0)4 74 66 69 50  
[export@carinamedical.com](mailto:export@carinamedical.com)