



Halogen lamp 50W

174 range

174 02

CARINA
Médical



CE
ISO 9001
ISO 13485

*Powerful and versatile, this lamp is ideal for general medicine
and minor surgery*

TECHNICAL DETAILS

- Halogen lamp 50 Watts, black colour
- Flexible arm, length 650mm
- Light intensity: 870 Lumens-3000° Kelvin (neutral white)
- CRI: 100/100
- Luminous power: 40000 Lux/ 50cm, illuminated diameter: 100mm
- Consumption: 50W
- Dichroic source 12V
- Bulb longevity: 5000h minimum
- Head diam. 11mm mounted on orientation joint
- Power cable length : 3m
- Epoxy covered base, white or black colour (RAL 9003 or RAL 9005)
- 5 rotating wheels
- Classe I medical device

OPTION

ISC (No contact switch)

You may also like :

Réf. **174 01**: Same as 174 02 but white colour

Réf. **176 02**: Same as 174 02 but wall-mounted

Réf. **176 01**: Same as 174 02 but white colour
and wall-mounted



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days
upon receipt of order
- ✓ Warranty: 5 years



CONTACT



Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com