



LED lamp 3W
178 range
178 02

CARINA
Medical



CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- LED lamp 3 Watts (beam 25°), black colour
- Flexible arm, length 650mm (also available in length 1000mm)
- Light intensity: 420 Lumens-4500° Kelvin (neutral white)
- CRI: 95/100 for an excellent colour rendering
- Luminous power: 7000 Lux/ 50cm, illuminated diameter: 250mm
- Low consumption: 3W
- Longevity: 50000h minimum (LED circuit, no bulb to replace)
- Head diam. 60mm mounted on orientation joint
- Power cable length : 3m
- Epoxy covered base, white colour or black colour (RAL 9003 or RAL 9005)
- 5 rotating wheels
- Classe I medical device

OPTION

ISC (no contact switch- Timer from 0 to 60 minutes)

**Ideal lamp for all types of care
and examinations in general
medicine, midwifery and
nursing**

You may also like :

Réf. **178 01**: Same as 178 02 but white colour

Réf. **178 04**: Same as 178 02 but wall-mounted

Réf. **178 03**: Same as 178 02 but white colour
and wall-mounted



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days
upon receipt of order
- ✓ Warranty: 5 years



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com