



Circular magnifying LED lamp 4W

179 range

179 01

CARINA
Medical



CE
ISO 9001
ISO 13485

The essential tool for dermatology or every other exams that requires meticulousness, finesse and precision

TECHNICAL DETAILS

- Circular magnifying LED lamp, white colour
- Flexible arm, length 860mm
- Light intensity: 300 Lumens-6400° Kelvin (neutral white)
- Luminous power: 2500 Lux/ 10cm
- Mineral glass lens diam. 120mm with magnification x2,25 (5 dioptries)
- Removable protective cap
- Class II isolation
- Delivered with clamp and weighted rolling base
- Power cable length : 1,5m
- Epoxy covering, white colour (RAL 9003)
- 5 rotating wheels



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 5 years



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com