



LED lamp 7W
178 range
188 02

CARINA
Médical



CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- LED lamp 7 Watts, black colour
- Flexible arm, length 650mm (also available in length 1140mm)
- Light intensity: 980 Lumens-4500° Kelvin (neutral white)
- CRI: 95/100 for an excellent colour rendering
- Luminous power: 60000 Lux/ 50cm, illuminated diameter: 180mm
- Low consumption: 7W
- Longevity: 50000h minimum (LED circuit, no bulb to replace)
- Power cable length : 3m
- Epoxy covered base, white colour or black colour (RAL 9003 or RAL 9005)
- 5 rotating wheels
- Classe I medical device

OPTION

ISC (no contact switch) or Dimmer

**Ideal lamp
for small surgery, general
medicine, dermatology and
podology**

You may also like :

Réf. **188 01**: Same as 188 02 but black colour



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 5 years



CONTACT



Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com