



LED lamp 17W

178 range

198 01

CARINA
Medical



CE
ISO 9001
ISO 13485

Very powerful treatment and intervention lamp, dedicated to small surgery, dermatology, intensive care or every other exams that requires homogeneous and precise lighting

TECHNICAL DETAILS

- High power LED lamp 17 Watts, white colour
- Flexible arm, length 700mm
- Light intensity: 1680 Lumens-4500° Kelvin (neutral white)
- CRI: 95/100 for an excellent colour rendering
- Luminous power: 120000 Lux/ 50cm, illuminated diameter: 200mm
- Consumption: 17W
- Longevity : 50000h minimum (LED circuit, no bulb to replace)
- Power cable length : 3m
- Epoxy covering, white colour (RAL 9003)
- 5 rotating wheels
- Classe I medical device

OPTION

ISC (no contact switch) and double intensity 50% / 100% switch



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 5 years



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com