



*Trolley*

**243 Range**

**243 08**

**CARINA**  
*Médical*



**CE**  
ISO 9001  
ISO 13485

## TECHNICAL DETAILS

- Trolley with 3 trays
- Dimensions: 400 x 400mm
- Square tube structure, epoxy coated
- With 2 side guards on the top of the first tray
- Equipped with 4 wheels diam 50mm

## DIMENSIONS

- Width: 40cm
- Length: 40cm
- Height: 80cm
- Load capacity: 50kg in static



### WHAT YOU NEED TO KNOW


- ✓ Deliver disassembled
- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 2 years



## BASE COLOURS




**RAL 9003**



**RAL 1013**



**RAL 1015**



**RAL 6019**



**RAL 6033**



**RAL 6000**



**RAL 5017**



**RAL 5002**



**RAL 5000**



**RAL 7035**



**RAL 7016**



**RAL 9005**




**RAL 1018**



**RAL 3020**



**RAL 3003**



**SILVER**



## MAINTENANCE GUIDE

### **DAILY MAINTENANCE**

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvavios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

### **WARRANTY**

Mechanical/ electrical part: 2 years



CONTACT



Parc d'entreprises Visionis  
278 rue de l'industrie  
01090 GUÉREINS- FRANCE

**COMMERCIAL SERVICE**

Tel: +33(0)4 74 66 69 49  
Fax: +33(0)4 74 66 69 50  
[export@carinamedical.com](mailto:export@carinamedical.com)