



Trolley
3070 Range
3071

CARINA
Médical



CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- Trolley with 2 trays
- Dimensions: 400 x 400mm
- Square tube structure, epoxy coated
- With one drawer
- Equipped with 4 wheels diam 50mm

DIMENSIONS

- Width: 40cm
- Length: 40cm
- Height: 80cm
- Load capacity: 50kg in static



WHAT YOU NEED TO KNOW


- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 2 years



BASE COLOURS




RAL 9003



RAL 1013



RAL 1015



RAL 6019



RAL 6033



RAL 6000



RAL 5017



RAL 5002



RAL 5000



RAL 7035



RAL 7016



RAL 9005




RAL 1018



RAL 3020



RAL 3003



SILVER



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvaniol pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part: 2 years



CONTACT



Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com