



Electric examination couch

326 range

326 08T

CARINA
Medical




CE
ISO 9001
ISO 13485

DESCRIPTION TECHNIQUE

- 3 sections examination couch, width 60cm
- Control : Pedal
- Backrest and legrest adjustment in every positions with gas jack
- Flat upholstery with good density 30kg/cm², M1 fire retardant class, thickness 8cm
- Equipped with paper roll holder for exam sheet, on the foot side
- Epoxy covering, light grey colour (RAL 7035) in standard

DIMENSIONS

- Width 60cm
- Adjustable height from 57 to 87cm
- Load capacity : 150kg in dynamic



You may also like :

Ref. **326 09T** : Same as 326 08T but width 70cm



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty : 5 years (1 year for the upholstery)

ACCESSORIES



ACC 997
4 retractable wheels diam. 50mm
for 326 electric height adjustment
range, 3 sections



ACC 990
Set of 4 wheels with brake
diam. 125mm



ACC 76B
Unit clamp to fix infusion stand
diam. 14mm for 326 electric
height adjustment range




ACC 76
Pair of clamps to fix infusion
stand diam. 14mm for 326
electric height adjustment range



ACC 76C
Unit clamp to fix infusion stand
diam. 16mm



ACC 32012
Stainless steel infusion stand
diam. 16mm, 2 plastic hooks,
without clamp



ACC 105
Pair of stainless steel stirrups
diam. 14mm, without clamp



ACC 964
Pair of legrests on knee system
in black polyurethane
diam. 14mm, without clamp



SLFEMIREPOSE
Cover for legrests (ACC964) in
upholstery, the unit




ACC 514-UNIVERSE
Plastic protection cover
for all couches



ACC 110POLY
Polyurethane blood test splint
without clamp, the unit



ACC 110
Blood test splint in colour
without clamp, the unit



ACC 32660
Half moon head rest width 60cm
fire retardant M1 class



ACC TETERGO
Tête ergonomique amovible



ACC 832
Face hole with plug cushion



ACC BARRILIT326B
Pair of retractable security side
rails for 326 with 3 sections

ACC 32614
Electric backrest for 326 electric
height adjustment range

ACC 32615
Electric legrest for 326 electric
height adjustment range

ZOOM
326 08T
with accessories

ACC 105
Pair of stainless steel stirrups
diam. 14mm



ACC 76
Pair of clamps diam. 14mm

UPHOLSTERY'S COLOURS

Fire retardant upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range
Be careful ! Urban colours **with extra value**

Upholstery : Classic range



Upholstery : Urban range







BASE COLOURS

The entire 326 range is offered in light grey (RAL 7035) as standard
Possibility to change the color without extra charge


Base colour




RAL 9003



RAL 1013



RAL 1015



RAL 6019



RAL 6033



RAL 6000




RAL 5017



RAL 5002




RAL 5000



RAL 7035



RAL 7016



RAL 9005




RAL 1018



RAL 3020



RAL 3003



Shiny

SILVER



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part : 5 years

Upholstery : 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL DEPARTMENT

Tél : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com