



Stool
40 range
4001ALU

CARINA
Medical



CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- Height adjustable stool with chromed gas column
- Manual control under the seat
- Fire retardant upholstery M1 class with visible seams, thickness 8cm
- Equipped with 5 wheels, diameter 50mm
- 5 branches aluminium base diameter 500mm

DIMENSIONS

- Round seat, diameter 35cm
- Adjustable height from 50 to 70cm
- Load capacity : 180kg

YOU MIGHT ALSO LIKE :

Réf. **4003ALU** : Same as 4001ALU but on fixed skids, adjustable height from 46 to 66cm



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 2 years (1 year for upholstery)

ACCESSORIES



ACC 10
Backrest for PRISCA,
range 120, 30 and 40



ACC 14
Foot control



ACC 15
Footrest



ACC 16
Chromed base,
diameter 38cm
(Reducing height of 1cm)



ACCROULETTE-02
5 single-wheels,
diameter 50mm



ACCROULETTE-03
5 double-wheels for hard
floors, diameter 50mm



ACCROULETTE-04
5 antistatic wheels, diameter
50mm



ACCROULETTE-05
5 auto-braking wheels,
in load or out of load

HH/HB
High or low height

UPHOLSTERY'S COLOURS

Fire retardant
upholstery M1 class

For the upholstery colours, you can choose between the Classic and Urban range

Upholstery : Classic range




Arno




Garonne



Potomac



Amour



Mississippi




St Laurent



Rhône



Missouri



Seine



Amazone

Upholstery : Urban Range



Abysse




Cendre



Origan




Géranium



Réglisse



Prune



Perle



Ficelle



Terre



Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate



MAINTENANCE GUIDE

DAILY MAINTENANCE

- Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material
- Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

- Mechanical / electrical part : 2 years
- Upholstery : 1 year



CONTACT



Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tél : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com