



*Fixed height sampling chair*

**512 range**

**512 02TB**

**CARINA**  
*Medical*




**CE**  
ISO 9001  
ISO 13485

## TECHNICAL DETAILS

- 3 sections sampling chair with fixed height 50cm
- Head rest adjustment in every positions with gas jack
- Legrest adjustment via a self locking mechanism system, the legrest maybe positioned vertically
- Trendelenburg position
- Flat upholstery with good density 35kg/cm<sup>2</sup>, M1 fire retardant class, thickness 8cm
- Equipped with 2 blood test splints, adjustable in every positions with serrated blocks system
- Epoxy covering, light grey colour (RAL 7035) in standard

## DIMENSIONS

- Width 55cm
- Height 50cm
- Load capacity : 130kg



## OTHER ELEMENTS

- All our chairs can be put in **lying position**
  - Low-height sampling chairs (50cm) are intended for **blood sampling**, higher chairs are dedicated to **gynecological sampling**.
- To give a better confort to patients and practitioners

### You may also like :

Ref. **512 02T** : Same as 512 02TB but height 67cm



## WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty : 2 years (1 year for the upholstery)

## ACCESSORIES




**ACC 537**  
Pair of clamps diam. 14mm for blood test splints and stirrups for all ranges except LUVIA



**ACC 535**  
Pair of double clamps diam. 14 and diam. 16mm to fix both armrest and infusion stand

**ACC 516**  
Option synchronized backrest and legrest



**ACC 105**  
Pair of stainless steel stirrups diam. 14mm without clamp



**ACC 964**  
Pair of legrests on knee system in black polyurethane, diam. 14 mm, without clamp



**SLFEMIREPOSE**  
Cover for legrests ACC964 in upholstery, the unit



**ACC 514-GOUT**  
Plastic protection for blood test splint, the unit




**ACC 514-ASSISE**  
Plastic protection for seatrest



**ACC 514-JAMBIÈRE**  
Plastic protection for legrest



**ACC 517-512**  
Plastic protection for backrest



**ACC 54005**  
Half moon head rest width 60cm fire retardant M1 class



**ACC 910-PLAT02**  
Stainless steel tray 340x340 with clamp


**ACC 910-PLAT04**  
Plastic tray 350x270 with clamp



**ACC 113U**  
Pair of handles without clamp with U shape



**ACC 113U-105**  
Pair of handles with clamps to fix both handles and stirrups



**ACC 125**  
Paper roll holder for all models



**ACC 32012**  
Stainless steel infusion stand diam. 16mm, 2 plastic hooks, without clamp



**ACC 32014**  
Stainless steel infusion stand diam. 16mm, 2 stainless steel hooks without clamp



**ACC 9505**  
Set of 4 wheels with brake diam. 50mm



**ACC 993-FIXE**  
4 retractable wheels diam. 50mm



**ACC TETERGO**  
Ergonomic headrest



**ACC 110**  
Blood test splint in colour without clamp, the unit

ZOOM  
PRODUCT

# ZQOM PRODUCT



*Photo of the model with plastic protections (option)*

# UPHOLSTERY'S COLOURS

Fire retardant upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range  
Be careful ! Urban colours **with extra value**

## *Upholstery : Classic range*




Arno



Garonne




Potomac



Amour




Mississippi



St Laurent



Rhône



Missouri



Seine



Amazone

## *Upholstery : Urban range*



Abysses



Cendre



Origan




Géranium




Réglisse




Prune




Perle



Ficelle




Terre



Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia




Tomate


# BASE COLOURS

The entire 512 range is offered in light grey (RAL 7035) as standard  
Possibility to change the color without extra charge


## *Base colour*




**RAL 9003**



**RAL 1013**



**RAL 1015**



**RAL 6019**



**RAL 6033**



**RAL 6000**




**RAL 5017**



**RAL 5002**




**RAL 5000**



**RAL 7035**



**RAL 7016**



**RAL 9005**




**RAL 1018**



**RAL 3020**



**RAL 3003**



*Shiny*

**SILVER**



# MAINTENANCE GUIDE

## **DAILY MAINTENANCE**

Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base : cleaning with soapy water only (no solvents) without abrasive materials

## **WARRANTY**

Mechanical/ electrical part : 2 years

Upholstery : 1 year



**CONTACT**

**CARINA**  
*Médical*

Parc d'entreprises Visionis  
278 rue de l'industrie  
01090 GUÉREINS- FRANCE

**COMMERCIAL DEPARTMENT**

Tél : +33(0)4 74 66 69 49  
Fax : +33(0)4 74 66 69 50  
[export@carinamedical.com](mailto:export@carinamedical.com)