



Hydraulic ENT chair

Range 512
51205T-ORL

CARINA
Medical




CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- 3 sections ENT chair with hydraulic height adjustment
- Rotation at 360° (possibility to block the rotation)
- Synchronised backrest and legrest adjustment, with gas jack
- Adjustable head rest specially developed for ENT exams
- Trendelenburg position
- Ergonomic upholstery with good density 30kg/cm², M1 fire retardant class, thickness 8cm (specially developed for an optimal comfort of the patient, with a lumbar support)
- Equipped with 2 removable armrests and a foot rest
- Epoxy covering, light grey color (RAL7035) in standard

DIMENSIONS

- Width 55cm
- Adjustable height from 54 to 74cm
- Load capacity: 130kg in dynamic



Adjustable head rest specially developed for ENT exams

ACCESSORIES



ACC514-ASSISE
ACC514-JAMBIÈRE
Plastic protection for seatrest or legrest



ACC 993
4 retractable wheels, diam.50mm



WHAT YOU NEED TO KNOW


- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 2 years (1 year for the upholstery)

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

*For the upholstery colours, you can choose between the Classic and Urban range
Be careful ! Urban colours with extra value*

Upholstery: Classic range



Arno




Garonne



Potomac




Amour



Mississippi



St Laurent




Rhône



Missouri



Seine




Amazone

Upholstery: Urban range



Abysses



Cendre



Origan




Géranium




Régilisse




Prune




Perle




Ficelle




Terre



Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate



**BASE
COLOURS**




RAL 9003



RAL 1013



RAL 1015



RAL 6019



RAL 6033



RAL 6000




RAL 5017



RAL 5002



RAL 5000



RAL 7035



RAL 7016



RAL 9005




RAL 1018



RAL 3020



RAL 3003



SILVER



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvaniol pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part: 2 years

Upholstery: 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com