



*Picture of the model with wheels option
(Ref. 52201-ORL+990)*

Fixed Height ENT chair

Range 522
52201-ORL

CARINA
Médical





CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- 3 sections ENT chair with fixed height
- Synchronised backrest and legrest adjustment, with gas jack
- Adjustable head rest specially developed for ENT exams
- Trendelenburg position
- Ergonomic upholstery with good density 30kg/cm³, M1 fire retardant class, thickness 8cm (specially developed for an optimal comfort of the patient, with a lumbar support)
- Equipped with 2 removable armrests
- Epoxy covering, light grey color (RAL7035) in standard

DIMENSIONS

- Width 55cm
- Height 50cm
- Load capacity: 130kg



Adjustable head rest
specially developed for
ENT exams

ACCESSORIES



ACC514-ASSISE
ACC514-JAMBIÈRE

Plastic protection for seatrest or legrest



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 2 years (1 year for the upholstery)

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

*For the upholstery colours, you can choose between the Classic and Urban range
Be careful ! Urban colours with extra value*

Upholstery: Classic range




Arno



Garonne



Potomac




Amour



Mississippi



St Laurent




Rhône



Missouri



Seine




Amazon

Upholstery: Urban range



Abyss



Cendre



Origan




Géranium




Régisse




Prune




Perle



Ficelle




Terre




Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia




Tomate


**BASE
COLOURS**



RAL 9003



RAL 1013



RAL 1015



RAL 6019



RAL 6033



RAL 6000



RAL 5017



RAL 5002



RAL 5000



RAL 7035



RAL 7016



RAL 9005




RAL 1018



RAL 3020



RAL 3003



SILVER



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part: 2 years

Upholstery: 1 year



CONTACT



Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com