



Fixed Height ENT chair

Range 522

52201-ORL + 990

CARINA
Medical



CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- 3 sections ENT chair with fixed height
- Synchronised backrest and legrest adjustment, with gas jack
- Adjustable head rest specially developed for ENT exams
- Trendelenburg position
- Ergonomic upholstery with good density 30kg/cm², M1 fire retardant class, thickness 8cm (specially developed for an optimal comfort of the patient, with a lumbar support)
- Equipped with 2 removable armrests
- Epoxy covering, light grey color (RAL7035) in standard
- On 4 wheels diam. 125 with brakes

DIMENSIONS

- Width 55cm
- Height 50cm
- Load capacity: 130kg



Adjustable head rest
specially developed for
ENT exams

ACCESSORIES



ACC514-ASSISE
ACC514-JAMBIÈRE
Plastic protection for seatrest or legrest

You may also like :

Réf. 522 01-ORL: Same as 522 01-ORL+990
but without wheels

Réf. 522 77T-ORL: Same as 522 01-ORL+990
but with chromed base, without wheels



WHAT YOU NEED TO KNOW


- ✓ Delivery time: 15 days
upon receipt of order
- ✓ Warranty: 2 years (1 year
for the upholstery)

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

*For the upholstery colours, you can choose between the Classic and Urban range
Be careful ! Urban colours with extra value*

Upholstery: Classic range



Arno




Garonne




Potomac




Amour



Mississippi



St Laurent



Rhône



Missouri



Seine




Amazone

Upholstery: Urban range




Abysses




Cendre




Origan




Géranium




Réglice




Prune




Perle




Ficelle




Terre




Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate



**BASE
COLOURS**




RAL 9003



RAL 1013



RAL 1015



RAL 6019



RAL 6033



RAL 6000




RAL 5017



RAL 5002



RAL 5000



RAL 7035



RAL 7016



RAL 9005




RAL 1018



RAL 3020



RAL 3003



SILVER



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvaniol pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part: 2 years

Upholstery: 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com