



Picture of the model with black ABS base instead of chromed base.

Polyurethane stool
60 range
6001CHR

CARINA
Medical



CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- Height adjustable stool with a chromed gas column
- Manual control under the seat
- Black polyurethane seat, thickness 6cm
- Equipped with 5 wheels, diameter 50mm
- 5 branches chromed base, diameter 550mm

DIMENSIONS

- Round seat diameter 35 cm
- Adjustable height from 49 to 69cm
- Load capacity : 130kg

You might also like :

Réf. **6011CHR** : Same as 6001CHR but on skids, adjustable height from 45 to 65cm



WHAT YOU NEED TO KNOW

- ✓ Delivery time : 15 days upon receipt of order
- ✓ Warranty : 2 years

ACCESSORIES



ACC 11
Black polyurethane backrest for
the 60 range



ACC 14
Foot control



ACC 16
Chromed base,
diameter 38cm
(Reducing height of 1cm)



ACC 15
Foot rest



ACCROULETTE-02
5 single-wheels,
diameter 50mm



ACCROULETTE-03
5 double-wheels for hard
floors, diameter 50mm



ACCROULETTE-04
5 antistatic wheels, diameter
50mm



ACCROULETTE-05
5 auto-braking wheels,
in load or out of load

HH/HB
High or low height



MAINTENANCE GUIDE

DAILY MAINTENANCE

- Base : Clean with soapy water only (solvents prohibited), without abrasive material

WARRANTY

- 2 years



CONTACT

CARINA
Medical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tél : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com