



Folding screen
881 range
881 58 13

CARINA
Medical





CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- Folding screens with 3 panels
- Panels mounted on hinges
- White PVC, fire retardant M1 class
- Chromed structure
- Equipped with 6 wheels diam. 50mm

DIMENSIONS

- Panels width 480x480x480mm
- Unfold length 1,44m
- Height 1,7m



WHAT YOU NEED TO KNOW

- ✓ Warranty: 3 years (1 year for the upholstery)

You may also like:

Réf. 881 58 12: Same as 881 58 13 but with 2 panels

Réf. 881 58 14: Same as 881 58 13 but with 4 panels



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

3 years

Upholstery: 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com