



Electric haemodialysis chair

945 range

945T

CARINA
Medical




CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- 3 sections haemodialysis chair, with electric height adjustment
- Multifunction hand remote control
- Adjustable backrest and legrest with electric jacks
- Tredelenburg position
- Upholstery with bultex foam, without visible seams, thickness 80mm
- Equipped with 2 blood test splints adjustable in every position by mechanical systems, infusion stand with 2 plastic hooks
- Welded base, all joint are mounted on bronze rings, height adjustment with electric jack
- Epoxy covered base, light grey colour (RAL 7035) in standard
- On 4 wheels with full brake diam.50mm

DIMENSIONS

- Width 55cm, width overall 85cm
- Adjustable height from 55 to 90cm
- Load capacity: 130kg in dynamic




WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 5 years (1 year for the upholstery)

ACCESSORIES



ACC 537
Pair of clamps diam.14mm
for blood test splints and stirrups



ACC 105
Pair of stainless steel stirrups
diam.14mm, without clamp



ACC 964
Pair of legrest on knee system
in black polyurethane,
diam.14mm, without clamp



SLFEMIREPOSE
Cover for legrests ACC964
in upholstery, the unit



ACC 514-GOUT
Plastic protection
For blood test splint, the unit




ACC 514-ASSISE
Protection plastic for seatrest



ACC 514-JAMBIÈRE
Plastic protection for legrest



ACC 517-512
Protection plastic for backrest



ACC 54005
Half moon headrest width 550mm
fire retardant M1 class



ACC 910-PLATO02
Stainless steel tray 340x340
with clamp

ACC 910-PLAT04

Plastic tray dimensions 350x270 with
clamps



ACC-113U
Pair of handles
With U shape



ACC 113U-105
Pair of handles with clamps to be
able to add a pair of stirrups on the
chair



ACC TETERGO
Ergonomic head rest




ACC 32012
Infusion stand diam.16mm with
2 polyamide hooks, without clamp



ACC 32014
Infusion stand diam.16mm,
2 stainless steel hooks, without clamp



ACC 110POLY
Polyurethane blood test splint
without clamp, the unit



ACC 125
Paper roll holder

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

*For the upholstery colours, you can choose between the Classic and Urban range
Be careful ! Urban colours with extra value*

Upholstery : Classic range




Arno



Garonne



Potomac



Amour



Mississippi



St Laurent




Rhône



Missouri



Seine




Amazone

Upholstery : Urban range



Aysse



Cendre



Origan




Géranium




Réglisse




Prune




Perle



Ficelle




Terre




Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia




Tomate




BASE COLOURS



RAL 9003



RAL 1013



RAL 1015



RAL 6019



RAL 6033



RAL 6000



RAL 5017



RAL 5002



RAL 5000



RAL 7035



RAL 7016



RAL 9005




RAL 1018



RAL 3020



RAL 3003



SILVER



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvantos pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part: 5 years

Upholstery: 1 year



CONTACT



Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com