



Photo of the model with all around controller and retractable wheels (Ref. 99906)

Electric examination couch

999 range

999 03

CARINA
Medical




CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- 2 sections examination couch, width 60cm
- Control : All around controller
- Backrest adjustment in every positions with gas jack
- Flat upholstery with good density 40kg/cm², M2 fire retardant class, thickness 7cm
- Frame composed of a 600N electric jack
- Equipped with paper roll holder for exam sheet, on the foot side
- Epoxy covering, silver colour
- On 4 skids

DIMENSIONS

- Width 60cm
- Adjustable height from 48 to 95cm
- Load capacity : 150kg in dynamic



You may also like :

Ref. **99901** : Same as 99903 but hand remote control

Ref. **99902** : Same as 99903 but pedal control

Ref. **99906** : Same as 99903 but on 4 retractable wheels

Ref. **99904** : Same as 99903 but hand remote control, on 4 retractable wheels

Ref. **99905** : Same as 99903 but pedal control, on 4 retractable wheels



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty : 2 years (1 year for the upholstery)

ACCESSORIES



ACC 538
Pair of clamps diam. 14mm
width 60cm




ACC 964
Pair of legrests on knee system
in black polyurethane
diam. 14mm, without clamp



ACC 76C
Unit clamp to fix infusion stand
diam. 16mm




ACC 110
Blood test splint in colour
without clamp, the unit



ACC 54001
Half moon head rest width
60cm fire retardant M1 class



ACC 514-UNIVERSE
Plastic protection cover
for all couches



ACC 105
Pair of stainless steel
stirrups diam. 14mm,
without clamp



SLFEMIREPOSE
Cover for legrests (ACC964)
in upholstery, the unit



ACC 32012
Stainless steel infusion stand
diam. 16mm, 2 plastic hooks,
without clamp



ACC 110POLY
Polyurethane blood test splint
without clamp, the unit



ACC TETERGO
Ergonomic head rest



ACC 832
Face hole with plug cushion

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M2 class

For the upholstery colours, you can choose between the Ecomax, Classic and Urban range
Be careful ! Classic and Urban colours **with extra value**

Upholstery : ECOMAX range*



Galapagos **Bahamas** **Blanc**



Tahiti **Malouines**

*Standard colors for ECOMAX, SIMPLEX and 999 range,
Classic and Urban colors available **on demand**

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range
Be careful ! Classic and Urban colours **with extra value**

Upholstery : Classic range




Arno




Garonne




Potomac



Amour




Mississippi



St Laurent



Rhône



Missouri



Seine




Amazone

Upholstery : Urban range



Abysses



Cendre



Origan




Géranium




Réglisse




Prune




Perle




Ficelle




Terre



Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part : 2 years

Upholstery : 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL DEPARTMENT

Tél : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com