



Folding massage table

WOOD range

ALUMINIUM FMA241

CARINA
Medical




CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- 1 section massage table
- Multi-adjustable headrest and face hole
- Vinyl upholstery, thickness 50mm
- With transport bag on wheels
- Aluminium frame

DIMENSIONS

- Width: 700mm
- Headrest length: 290mm
- Adjustable height from 68 to 92cm
- Load capacity: 200kg
- Table weight: 12kg (**very light table, easily transportable**)



Upholstery colors




Beige

Blue

Chocolate

Black



WHAT YOU NEED TO KNOW



Warranty: 1 year

ACCESSORIES

CUSHIONS



COUSSIN-001

Rectangular cushion with HR foam
dim. 440x200x75mm




COUSSIN-002

Half cylindric cushion HR foam,
dim 450x120mm, radius 120mm



COUSSIN-003

Cylindric cushion with HR foam
dim. 500x160mm



COUSSIN-004

Triangular cushion with HR foam
dim. 600x500x300mm



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvaniol pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com