



Physiotherapy table

BOBATH range

BOBATH2

CARINA
Medical




CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- 1 section physiotherapy table
- Stretched upholstery, M1 fire retardant class, thickness 4 cm
- Epoxy covered base, light grey color (RAL7035) in standard

DIMENSIONS

- Width: 1000mm
- Length: 2000mm
- Height: 500mm
- Load capacity: 150kg



You may also like :

Réf. BOBATH : Same as BOBATH2 but width 2000mm




WHAT YOU NEED TO KNOW


- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 2 years (1 year for the upholstery)

ACCESSORIES

CUSHIONS




COUSSIN-001
Rectangular cushion with HR foam
dim. 440x200x75mm



COUSSIN-002
Half cylindric cushion HR foam,
dim 450x120mm, radius 120mm



COUSSIN-003
Cylindric cushion with HR foam
dim. 500x160mm




COUSSIN-004
Triangular cushion with HR foam
dim. 600x500x300mm

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

*For the upholstery colours, you can choose between the Classic and Urban range
Be careful ! Urban colours with extra value*

Upholstery: Classic range



Arno




Garonne




Potomac




Amour



Mississippi



St Laurent




Rhône



Missouri



Seine




Amazone

Upholstery: Urban range



Abysses



Cendre



Origan




Géranium




Régilisse




Prune




Perle




Ficelle




Terre




Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate



BASE COLOURS

*All the BOBATH range have a base grey in standard
Possibility to change without extra value*



RAL 9003



RAL 1013



RAL 1015



RAL 6019



RAL 6033



RAL 6000




RAL 5017



RAL 5002



RAL 5000



RAL 7035



RAL 7016



RAL 9005




RAL 1018



RAL 3020



RAL 3003



SILVER



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvaniol pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part: 2 years

Upholstery: 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com