



Picture of the model with pedal

Aesthetic chairs

INJEXIA range

INJEXIA05

CARINA
Medical



CE


ISO 9001
ISO 13485

TECHNICAL DETAILS

- 2 sections aesthetic chair, with electric height adjustment
- Control: Hand remote
- Backrest adjustment with electric jack
- Equipped with 2 armrests synchronized with backrest and an adjustable headrest
- Trendelenburg position
- Design upholstery with a good density, viscoelastic foam, multi-layer coating guaranteeing longevity and resistance to tearing and friction, thickness 80mm
 - Anti-fouling treatment
 - Fungicide treatment
 - Antimicrobial treatment
- Epoxy covering, sand color in standard
- With 2 small wheels on the head side and 2 adjustable skids on the foot side

DIMENSIONS

- Width 60cm
- Adjustable height from 48 to 88cm
- Load capacity: 200kg in dynamic



You may also like :

Réf. **INJEXIA06**: Same as INJEXIA05 but multifunctions pedal control



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 5 years (1 year for the upholstery)

ACCESSORIES




ACCINJEXIA01
Tray support with tray



ACCINJEXIA04
Leg cushion to transform into table



ACCINJEXIA02
Stems for lamp's support
without lamp



ACC CLAMEDIAM16
Adaptation clamp in aluminium
for fixation on rail dim25x10, the unit
diam. 14mm, diam. 16mm




No inclination of the legrest.
Possibility to make the chair horizontal
by adding a cushion **ACCINJEXIA04**

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range

Upholstery: Classic range



Arno




Garonne



Potomac




Amour




Mississippi




St Laurent



Rhône



Missouri



Seine




Amazone

Upholstery: Urban range




Abysses




Cendre




Origan




Géranium




Réglice




Prune




Perle




Ficelle




Terre



Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate




**BASE
COLOURS**




RAL 9003



RAL 1013



RAL 1015



RAL 6019



RAL 6033



RAL 6000




RAL 5017



RAL 5002



RAL 5000



RAL 7035



RAL 7016



RAL 9005




RAL 1018



RAL 3020



RAL 3003



SILVER



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvaniol pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

Mechanical/ electrical part: 5 years

Upholstery: 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com