



ENT CHAIR  
**OREA RANGE**  
**OREA-HVH1-SERENITY**

**CARINA**  
*Medical*



**CE**  
ISO 9001  
ISO 13485




## TECHNICAL DETAILS

- ENT chair with hydraulic height adjustment, 3 sections, on chromed base
- Equipped with an height adjustable headrest, also adjustable in inclination. Perfectly designed for ENT exams
- Chair is rotative at 360°
- Backrest and legrest/foot rest adjustments are synchronised with a gas jack
- Supplied with 2 adjustable arm rests
- Trendelenburg position
- SERENITY upholstery: Flat without visible seams- fire retardant, M1 class
- Chromed base

## DIMENSIONS

- Width 55 cm in standard or width 65 cm in option
- Height adjustment from 48 to 88cm
- Dynamic load capacity :  
150 Kg in standard configuration  
200 Kg if choose of electric backrest and legrest options



Adjustable headrest in height and inclination  
Specially designed for ENT exams



### You may also like:

Ref. **OREA-HF** : Fixed height

Ref. **OREA-HVE** : Electric height



## WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days from date of order
- ✓ Warranty : 2 years (1 year for the upholstery)

# ACCESSORIES



**SERENITY**  
Flat  
Upholstery



**CONFORT**  
Shaped  
upholstery



**PREMIUM**  
Visco elastic  
foam  
upholstery

## **ACC 514-JAMBIERE**

Plastic protection of  
the legrest - for chairs  
width 55cm

## **ACC 514-ASSISE**

Plastic protection of the  
the seat rest -for chairs  
width 55cm



## **ACC 993**

Set of 4 removable wheels  
diam. 50mm

# UPHOLSTERY'S COLOURS

Fire retardant,  
M1 class

## *Upholstery : Classic range*




## *Upholstery : Urban range*



## *Upholstery : GINKGO range*





















**BASE  
COLOURS**

All the range is suggested with a light grey colour RAL 7035.  
Possibility to change this colour without extra costs

*Epoxy Painting colours*

				
<b>RAL 9003</b>	<b>RAL 1013</b>	<b>RAL 1015</b>	<b>RAL 6019</b>	<b>RAL 6033</b>
				
<b>RAL 6000</b>	<b>RAL 5017</b>	<b>RAL 5002</b>	<b>RAL 5000</b>	<b>RAL 7035</b>
				
<b>RAL 7016</b>	<b>RAL 9005</b>	<b>RAL 1018</b>	<b>RAL 3020</b>	<b>RAL 3003</b>
				
	<b>SILVER</b>			



## MAINTENANCE GUIDE

### DAILY MAINTENANCE

Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base : cleaning with soapy water only (no solvents) without abrasive materials

### WARRANTY

- Structure : 2 years
- Upholstery : 1 year



CONTACT

CARINA  
*Médical*

Parc d'entreprises Visionis  
278 rue de l'industrie  
01090 GUÉREINS - FRANCE  
**COMMERCIAL SERVICE**  
Tél : +33(0)4 74 66 69 49

Fax : +33(0)4 74 66 69 50  
[export@carinamedical.com](mailto:export@carinamedical.com)