



Cardiology couch

**OVALIA V2 range**

**OV2-98-125**

**CARINA**  
*Medical*




**CE**  
ISO 9001  
ISO 13485

## TECHNICAL DETAILS

- 2 sections cardiology couch, width 80cm
- Control : Multifonctionnal hand remote
- Head rest adjustment in every positions with electric jack
- Trendelenburg and anti tredelenburg position in standard
- Flat upholstery with good density 24kg/m<sup>2</sup>, M1 fire retardant class, thickness 5cm
- Frame composed of 2 electronic columns, giving the couch high stability
- Equipped with cardiology removable section to make heart echography or other exams and a paper roll holder for exam sheet rollet on the head side
- Shapped base for ideal cleaning, base colour light grey (RAL 7035)
- On 4 wheels diam. 125mm with one brake on each wheel

## DIMENSIONS

- Width 80cm
- Adjustable height from 48 to 88cm
- Load capacity : 200kg in dynamic



Cardiology section removable to facilitate heart echography and other chest examinations

### You may also like :

Ref. **OV2-95-125** : Same as OV2-98-125 but width 70cm

Ref. **OV2-99-125** : Same as OV2-98-125 but multifonctionnal pedal control

Ref. **OV2-96-125** : Same as OV2-98-125 but width 70cm and multifonctionnal pedal control



## WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty : 5 years (1 year for the upholstery)

# ACCESSORIES



**RL-150**  
4 wheels diam. 150mm  
with central brake system



**ACC 76-OVALIA-80**  
Pair of clamps  
for couch width 80cm




**ACC 110POLY**  
Polyurethane blood test splint  
without clamp, the unit



**ACC 110**  
Blood test splint in colour  
without clamp, the unit



**ACC 964**  
Pair of legrests on knee system  
in black polyurethane,  
diam. 14 mm, without clamp



**ACC 105**  
Pair of stainless steel stirrups  
diam. 14mm, without clamp



**SLFEMIREPOSE**  
Cover for legrests ACC964  
in upholstery, the unit



**ACC 514-UNIVERSE**  
Plastic protection cover  
for all couches




**ACC 32012**  
Stainless steel infusion stand  
diam. 16mm, 2 plastic hooks,  
without clamp



**ACC 76C**  
Unit clamp to fix infusion stand  
diam. 16mm



**ACC 832**  
Face hole with plug cushion



**ACC 54003**  
Half moon headrest width 80cm  
fire retardant M1 class



**ACC OVALIA06**  
Hand remote for OVALIA +  
connecting cable for double  
control)



**ACC OVALIA05**  
Pedal for OVALIA + connecting  
cable for double control



**ACC OVALIA02**  
Pair of stainless steel rails 25x10  
length 1500mm for OVALIA V2  
2 sections



**ACC BATTERIE01**  
Battery, capacity 20 cycles  
(1 cycle = 1x height's up/ down +  
1x backrest's up/ down)




**ACC BARRILITOV08**  
Pair of retractable security side  
rails for couch width 80cm



**ACC TETERGO**  
Ergonomic headrest



**ACC-CAPOT**  
Base colour in a different colour  
than the standers light grey  
available in different colours



**ACC FEMINA01**  
Coloscope holder

# UPHOLSTERY'S COLOURS

Fire retardant  
upholstery, M1 class

For the upholstery colours, you can choose between the Classic and Urban range  
Be careful ! Urban colours **with extra value**

## *Upholstery : Classic range*




Arno



Garonne




Potomac



Amour




Mississippi



St Laurent




Rhône



Missouri



Seine



Amazone

## *Upholstery : Urban range*



Abysses



Cendre



Origan




Géranium




Régilisse




Prune




Perle




Ficelle




Terre




Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia




Tomate


# BASE COLOURS

The entire OVALIA V2 range has a light grey (RAL 7035) cover as standard  
Possibility to change the colour of the cover **with an extra charge (option ACC CAPOT)**


## *Base colour*




RAL 9003



RAL 1013



RAL 1015



RAL 6019



RAL 6033



RAL 6000




RAL 5017



RAL 5002




RAL 5000



RAL 7035



RAL 7016



RAL 9005




RAL 1018



RAL 3020



RAL 3003



*Shiny*

SILVER



# MAINTENANCE GUIDE

## **DAILY MAINTENANCE**

Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base : cleaning with soapy water only (no solvents) without abrasive materials

## **WARRANTY**

Mechanical/ electrical part : 5 years

Upholstery : 1 year



**CONTACT**

**CARINA**  
*Médical*

Parc d'entreprises Visionis  
278 rue de l'industrie  
01090 GUÉREINS- FRANCE

**COMMERCIAL DEPARTMENT**

Tél : +33(0)4 74 66 69 49  
Fax : +33(0)4 74 66 69 50  
[export@carinamedical.com](mailto:export@carinamedical.com)