



Stool
PRISCA range
PRISCA/ONE

CARINA
Medical



CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- Height adjustable stool with black column
- Manual control under the seat
- Fire retardant upholstery M2 class, thickness 5cm
- Equipped with 5 wheels, diameter 50mm
- 5 branches black ABS base, diameter 550mm

DIMENSIONS

- Round seat, diameter 35cm
- Adjustable height from 50 to 70cm
- Load capacity : 130kg



You might also like :

Ref. PRISCA/TWO : Same as PRISCA/ONE but with seatrest thickness 8cm, visible seams and height adjustment from 53 to 73cm



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days upon receipt of order
- ✓ Warranty: 2 years (1 year for upholstery)

UPHOLSTERY'S COLOURS

Fire retardant
upholstery, M2 class

For upholstery colours, you have the choice between the ECOMAX, Classic and Urban range.
Attention ! Classic and Urban colours with **added value**


Upholstery : ECOMAX range*



Galapagos



Bahamas



Blanc



Tahiti



Malouines

*Standard colours for ECOMAX and SIMPLEX range,
Classic and Urban colours available **on request**


UPHOLSTERY'S COLOURS

Fire retardant upholstery, M1 class

Upholstery : Classic range



Arno




Garonne




Potomac



Amour




Mississippi



St Laurent



Rhône



Missouri



Seine




Amazone

Upholstery : Urban Range



Abysses




Cendre



Origan




Géranium




Réglisse




Prune



Perle




Ficelle



Terre




Béton



Chataigne



Chocolat



Pensée



Pétrole



Turquoise



Granny



Absinthe



Ananas



Orange



Fuschia



Tomate



MAINTENANCE GUIDE

DAILY MAINTENANCE

- Upholstery : Clean with soapy water only (solvents prohibited), can be disinfected with products such as Phagosept Spray Salvanios pH10 at 0.75% and Anios Surfaces Hautes without abrasive material
- Base : Clean with soapy water only (solvents prohibited), without abrasive material
matériau abrasif

WARRANTY

- Base : 2 years
- Upholstery : 1 year



CONTACT

CARINA
Medical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tél : +33(0)4 74 66 69 49
Fax : +33(0)4 74 66 69 50
export@carinamedical.com