



Picture of the product with CONFORT type of upholstery

Blood sampling chair

SANA range

SANA-HF

CARINA
Medical




TECHNICAL DETAILS

- 3 sections blood sampling chair with fixed height
- Backrest adjustment with gas jack
- Legrest adjustment with quick, secured and easily adjustable mechanical system
- Trendelenburg position
- Large range of available upholsteries — SERENITY, CONFORT and PREMIUM — All of those certified fire retardant, M1 class
- Equipped with a pair of adjustable arm rests with indexable finger
- Equipped of 4 adjustable skids to optimize the chair's stability
- Epoxy base, light grey colour RAL 7035 in standard

DIMENSIONS

- Width 55 cm or 65 cm (in option)
- Height of 48 cm
- Dynamic load capacity :
150 Kg if bought with gas jack backrest
200 Kg if bought with electric backrest and legrest options



You may also like:

Réf. **SANA-HVE** Same as SANA-HF but with electric height adjustment

Réf. **SANA-HVH** Same as SANA-HF but with hydraulic height adjustment



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days from date of order
- ✓ Warranty : 5 years
1 year for the upholstery)

ACCESSORIES



SERENITY
Flat upholstery



CONFORT
Shaped upholstery



PREMIUM
Visco elastic foam upholstery



RL-150
4 directionnal wheels diam 150 mm , with central brake system To move the couch with a patient on



RL-125
4 directionnal wheels diam 150 mm , with central brake system To move the couch with a patient on



RL-50
4 wheels diam 50 mm with individuals breaks To move the couch with nobody on



ACC444-DELEC-1
Electrically adjustable backrest, for chair with fixed height



ACC444-JELEC-1
Electrically adjustable legrest, for chair with fixed height



ACC444-PEDALE
Pedal + Y cable to connect both pedal and hand remote - To be order with electric backrest and legrest options



ACC 444-TCDE
Hand remote + Y cable to connect both hand remote and pedal - To be order with electric backrest and legrest options



ACC 110-ROTULE1
Arm rest for blood test, with upholstery, on ball joint adjustment system, the unit , without clamp



ACC 110
Arm rest for blood test, with upholstery, the unit, without clamp



ACC 110-POLY
Arm rest for blood test in black polyurethane, the unit, without clamp



ACC TECH114
Pair of arm rests with one hand adjustment, without clamps



ACC 537
Pair of clamps diam 14mm for arm rests and stirrups, for chairs width 55 cm



ACCACCOU01
Pair of removable arm rests in polyuréthane



ACC 537-650
Pair of clamps diam 14mm for arm rests and stirrups, for chairs width 55 cm

ACCESSORIES



ACC 105

Pair of standard stainless steel stirrups, diam 14mm, without clamp



SLFEMIREPOSE

Cover for leg rest ref ACC964, Upholstery colour to be chosen, the unit



ACC BAVETTE01

Cover upholstery part between seat rest and leg rest



ACCBACGYNECO

Removable stainless steel tray for gynaecology



ACCMUL-320INO1

Infusion stand with plastic hooks, diam 16mm, without clamp

ACC 514-ASSISE

Plastic protection for seatrest, the unit

ACC 514-JAMBIERE

Plastic protection for legrest of chair width 55cm



ACC TETERGO

Têteière ergonomique amovible



ACC 125

Porte rouleau



ACC 964

Pair of leg rests, on ball-joint adjustment system, in black polyurethane, diam 14mm, without clamp



ACC080

Stainless steel tray

ACC113U

Pair of handles to help to sit and get up, without clamps, ready for use

ACC 514-GOUT

Plastic protection for arm rest, the unit

ACC 517-512

Plastic protection for backrest, the unit



ACC 540 05

Half moon headrest, width 550 mm, fire retardant M1 class



ACCCOUSSIN-01

Pillow cushion dim 500 x 400, fire retardant, M1 class

UPHOLSTERY COLOURS

Fire retardant,
M1 class

Upholstery : Classic range




Upholstery : Urban range



Upholstery : GINKGO range










**BASE
COLOURS**

All the range is suggested with a light grey colour RAL 7035. Possibility to change this colour without extra costs

Epoxy Painting colours

				
RAL 9003	RAL 1013	RAL 1015	RAL 6019	RAL 6033
				
RAL 6000	RAL 5017	RAL 5002	RAL 5000	RAL 7035
				
RAL 7016	RAL 9005	RAL 1018	RAL 3020	RAL 3003
				
SILVER				




MAINTENANCE GUIDE

DAILY MAINTENANCE

- Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvaniol pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material
- Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

- Structure: 5 years
- Upholstery : 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie 01090
GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tél: +33(0)4 74 66 69 49

Fax: +33(0)4 74 66 69 50
commercial@carinamedical.com