



*Picture of the model with CONFORT upholstery*

*Blood sampling chair*  
**SANA range**  
**SANA-HVE**

**CARINA**  
*Medical*




## TECHNICAL DETAILS

- 3 sections blood sampling chair with electric height adjustment, with pedal or hand remote
- Backrest adjustment with gas jack
- Legrest adjustment with quick, secured and easily adjustable mechanical system  
Trendelenburg position
- Large range of available upholsteries — SERENITY, CONFORT and PREMIUM — All of those certified fire retardant, M1 class
- Equipped with a pair of adjustable arm rests with indexable finger
- Equipped of 4 adjustable skids to optimize the chair's stability
- Epoxy base, light grey colour RAL 7035 in standard

## DIMENSIONS

- Width 55 cm or 65 cm
- Height adjustment from 48 cm to 88 cm
- Dynamic load capacity :  
150 Kg if bought with gas jack backrest  
200 Kg if bought with electric backrest and legrest options



### You may also like:

Réf. **SANA-HF**: Same as **SANA-HVE** but with fixed height

Réf. **SANA-HVH**: Same as **SANA-HVE** but with hydraulic height adjustment



## WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days from date of order
- ✓ Warranty : 5 years (1 year for the upholstery)

# ACCESSORIES



**SERENITY**  
Flat upholstery



**CONFORT**  
Shaped upholstery



**PREMIUM**  
Visco elastic foam upholstery



**RL-150**  
4 directionnal wheels diam 150 mm , with central brake system  
To move the couch with a patient on



**RL-125**  
4 wheels diam 125 mm with individuals breaks  
To move the couch with a patient on ( on small distances)



**RL-50**  
4 wheels diam 50 mm with individuals breaks  
To move the couch with nobody on



**ACC444-DELEC-2**  
Electrically adjustable backrest, for chair with electric height adjustment



**ACC444-JELEC-2**  
Electrically adjustable legrest, for chair with electric height adjustment



**ACC 444-TCDE**  
Hand remote + Y cable to connect both hand remote and pedal - To be order with electric backrest and legrest options



**ACC 110**  
Arm rest for blood test, with upholstery, the unit, without clamp



**ACC 444-PEDALE**  
Pedal + Y cable to connect both pedal and hand remote - To be order with electric backrest and legrest options



**ACC TECH114**  
Pair of arm rests with one hand adjustment, without clamps



**ACC 110-ROTULE1**  
Arm rest for blood test, with upholstery, on ball joint adjustment system, the unit , without clamp



**ACCACCOU01**  
Pair of removable arm rests in polyuréthane



**ACC 110POLY**  
Arm rest for blood test, in black polyurethane, the unit

## ACCESSORIES



### ACCCOUSSIN-01

Pillow cushion dim 500 x 400, fire retardant, M1 class



### ACC TETERGO

Ergonomic and removable headrest



### ACC 540 05

Half moon headrest, width 550 mm, fire retardant M1 class



### ACC 537

Pair of clamps diam 14mm for arm rests and stirrups, for chairs width 55 cm



### ACC 537-650

Pair of clamps diam 14mm for arm rests and stirrups, for chairs width 65 cm



### ACC 964

Pair of leg rests, on ball-joint adjustment system, in black polyurethane, diam 14mm, without clamp



### SLFEMIREPOSE

Cover for leg rest ref ACC964. Upholstery colour to be chosen, the unit



### ACC080

Stainless steel tray



### ACC BAVETTE01

Cover upholstery part between seat rest and leg rest

### ACC113U

Pair of handles to help to sit and get up, without clamps, ready for use



### ACCMUL-320INO1

Infusion stand with plastic hooks, diam 16mm, without clamp

### ACC 514-GOUT

Plastic protection for arm rest, the unit

### ACC 514-ASSISE

Plastic protection for seatrest, the unit

### ACC 517-512

Plastic protection for backrest, the unit

### ACC 514-JAMBIERE

Plastic protection for leg rest, the unit

# UPHOLSTERY COLOURS

Fire retardant,  
M1 class

## *Upholstery : Classic range*




## *Upholstery : Urban range*



## *Upholstery : GINKGO range*





**BASE  
COLOURS**

*All the range is suggested with a light grey colour RAL 7035. Possibility to change this colour without extra costs*

*Epoxy Painting colours*




**RAL 9003**



**RAL 1013**




**RAL 1015**




**RAL 6019**



**RAL 6033**



**RAL 6000**



**RAL 5017**



**RAL 5002**




**RAL 5000**



**RAL 7035**



**RAL 7016**



**RAL 9005**



**RAL 1018**



**RAL 3020**



**RAL 3003**



**SILVER**


# MAINTENANCE GUIDE

## DAILY MAINTENANCE

- Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material
- Base : cleaning with soapy water only (no solvents) without abrasive materials

## WARRANTY

- Structure : 5 years
- Upholstery : 1 year



CONTACT

CARINA  
*Médical*

Parc d'entreprises Visionis  
278 rue de l'industrie 01090  
GUÉREINS- FRANCE

**COMMERCIAL SERVICE**

Tél: +33(0)4 74 66 69 49

Fax: +33(0)4 74 66 69 50  
[export@carinamedical.com](mailto:export@carinamedical.com)