



Picture without CONFORT type of upholstery

Blood sampling chair
SANA Range
SANA-HVH2

CARINA
Medical




TECHNICAL DETAILS

- 3 sections blood sampling chair with hydraulic height adjustment, swivelling at 360°
- Backrest adjustment with gas jack
- Legrest adjustment with quick, secured and easily adjustable mechanical system
- Trendelenburg position
- Large range of available upholsteries — SERENITY, CONFORT and PREMIUM — All of those certified fire retardant, M1 class
- Equipped with a pair of adjustable arm rests with indexable finger
- With 4 fixed skids for chair's stability
- Epoxy base, light grey colour RAL 7035 in standard

DIMENSIONS

- Width 55 cm or 65 cm (in option)
- Height adjustment from 48 cm to 68 cm
- Dynamic load capacity :
150 Kg if bought with gas jack backrest
200 Kg if bought with electric backrest and legrest options



You may also like:

Réf. **SANA-HF**: Same as SANA-HVE but with fixed height

Réf. **SANA-HVE**: Same as SANA-HVH but with electric height adjustment



WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days from date of order
- ✓ Warranty : 5 years
- ✓ 1 year for the upholstery)

ACCESSORIES



SERENITY
Flat
upholstery



CONFORT
Shaped
upholstery



PREMIUM
Visco elastic
foam
upholstery



RL-150
4 roulettes diamètre 150 mm
directionnelles à freinage
centralisé Déplacement du
fauteuil avec patient Hauteur de
48 cm



RL-125
4 Roulettes diamètre 125
mm à freins indépendants
Déplacement du fauteuil
avec patient Hauteur de
48 cm



RL-50
4 roulettes diamètre 50 mm à
freins indépendants
Déplacement du fauteuil seul
Hauteur de 48 cm



ACC 110-ROTULE1
Arm rest for blood test,
with upholstery, on ball
joint adjustment system,
the unit, without clamp



ACC TECH114
Pair of arm rests with one
hand adjustment, without
clamps



ACC 110 POLY
Arm rest for blood test, in
black polyurethane, the
unit



ACCACCOU01
Pair of removable
arm rests in
polyuréthane



ACCCOUSSIN-01
Pillow cushion dim 500 x
400, fire retardant, M1
class



ACC 540 05
Half moon headrest,
width 550 mm, fire
retardant M1 class



ACC TETERGO
Ergonomic and
removable headrest



ACC 537-650
Pair of clamps diam
14mm for arm rests
and stirrups, for chairs
width 65 cm



ACC 537
Pair of clamps diam
14mm for arm rests
and stirrups, for chairs
width 55 cm

ACCESSOIRES



ACC 105

Pair of standard stainless steel stirrups, diam 14mm, without clamp ACC 105



SLFEMIREPOSE

Cover for leg rest ref ACC964, Upholstery colour to be choosed, the unit



ACC BAVETTE01

Cover upholstery part between seat rest and leg rest



ACCBACGYNECO

Removable stainless steel tray for gynaecology



ACCMUL-320INO1

Infusion stand with plastic hooks, diam 16mm, without clamp

ACC 514-ASSISE

Plastic protection for seatrest, the unit

ACC 517-512

Plastic protection for backrest, the unit



ACC 125

Paper roll holder



ACC 964

Pair of leg rests, on ball-joint adjustment system, in black polyurethane, diam 14mm, without clamp



ACC080

Stainless steel tray

ACC113U

Pair of handles to help to sit and get up, without clamps, ready for use

ACC 514-GOUT

Plastic protection for arm rest, the unit

ACC 514-JAMBIERE

Plastic protection for leg rest, the unit

UPHOLSTERY COLOURS

Fire retardant,
M1 class

Upholstery : Classic range




Upholstery : Urban range



Upholstery : GINKGO range









**BASE
COLOURS**

All the range is suggested with a light grey colour RAL 7035. Possibility to change this colour without extra costs

Epoxy Painting colours

				
RAL 9003	RAL 1013	RAL 1015	RAL 6019	RAL 6033
				
RAL 6000	RAL 5017	RAL 5002	RAL 5000	RAL 7035
				
RAL 7016	RAL 9005	RAL 1018	RAL 3020	RAL 3003
				
SILVER				


MAINTENANCE GUIDE

DAILY MAINTENANCE

- Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvaniol pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material
- Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

- Structure: 5 years
- Upholstery : 1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie 01090
GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tél: +33(0)4 74 66 69 49

Fax: +33(0)4 74 66 69 50
export@carinamedical.com