



Folding massage chair

WOOD range

WOODCHAISE FMA005

CARINA
Medical




CE
ISO 9001
ISO 13485

TECHNICAL DETAILS

- Folding massage chair
- Vinyl upholstery, thickness 50mm
- With adjustable chest cushion, seat, folding headrest and armrest
- With transport bag
- Folding articulated aluminium frame

DIMENSIONS

- Minimum height 500 tiltable
- Load capacity: 130kg
- Chair weight: 7kg



Upholstery colors



Beige




Blue



Chocolate



Black



WHAT YOU NEED TO KNOW



Warranty: 1 year



MAINTENANCE GUIDE

DAILY MAINTENANCE

Upholstery: Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvaniol pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base: cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

1 year



CONTACT

CARINA
Médical

Parc d'entreprises Visionis
278 rue de l'industrie
01090 GUÉREINS- FRANCE

COMMERCIAL SERVICE

Tel: +33(0)4 74 66 69 49
Fax: +33(0)4 74 66 69 50
export@carinamedical.com