

Electric examination couch

XENA range

XENA-2P-PREMIUM

CARINA Médical

MADE IN FRANCE

CE

ISO 9001
ISO 13485


BRONZE | Top 35%
ecovadis
Sustainability Rating
AUG 2025


TECHNICAL DETAILS

- Examination couch with eletrical height adjustment, 2 sections
- Control : Hand remote, pedal or pedal with memorization functions
- Backrest adjustment in every position by gas jack
- Premium upholostery, fire retardant M1 class, thickness 8 cm : very comfortable
- Equipped with paper roll holder on the foot side
- Base especially designed for optimal cleaning, colour light grey RAL 7035

DIMENSIONS

- Width 60 or 70 cm
- Height adjustment from 48 to 88cm
- Load capacity :
  150 kg if buyed with gas jack backrest
  200 kg if buyed with electric backrest option





You may also like:

Réf. **XENA-3P-PREMIUM :** Same as XENA-2P-PREMIUM but in 3 sections

WHAT YOU NEED TO KNOW

- ✓ Delivery time: 15 days from date of order
- ✓ Warranty : 5 years (1 year for the upholstery)

ACCESSORIES



RL-150
4 directionnal wheels diam 150 mm with central brake system
To move the couch with a patient on



RL-50
4 wheels diam 50 mm with individuals breaks
To move the couch with nobody on



ACC 444-TCDE
Hand remote + Y cable to connect both hand remote and pedal - To be order with electric backrest option



ACC 76
Pair of clamps to put arm rest or stirrups diam 14mm



ACC 110POLY
Arm rest for blood test in black polyurethane, the unit



ACC 964
Pair of leg rests, on ball-joint adjustment system, in black polyurethane, diam 14mm, without clamp



SLFEMIREPOSE
Cover for leg rest ref ACC964
Upholstery colour to be choosed, the unit



ACC 32012
Infusion stand with 2 plastic hooks, diam 16mm, the unit




ACC 832
Face hole with cushion plug



ACC BARRILITOV02
Pair of removable security side rails for exam couches width 60 or 70cm



RL-125
4 wheels diam 125 mm with individuals breaks
To move the couch with a patient on (on small distances)




ACC444-DELEC-1
Electrically adjustable backrest, on couch with electric height adjustment



ACC 444-PEDALE
Pedal + Y cable to connect both pedal and hand remote - To be order with electric backrest option




ACC 110
Arm rest for blood test, with upholstery, the unit



ACC 105
Pair of standard stainless steel stirrups, diam 14mm, without clamp



ACC 514-UNIVERSE
Transparent universal plastic protection for leg/foot section



ACC 54002
Half moon headrest, width 600mm, fire retardant M1 class



ACC TETERGO
Ergonomic and removable headrest

UPHOLSTERY COLOURS

Fire retardant, M1 class

Upholstery : Classic range



Upholstery : Urban range



Upholstery : GINKGO range



BASE COLOURS

All the range is suggested with a light grey colour RAL 7035. Possibility to change this colour with extra costs

Epoxy Painting colours



DAILY MAINTENANCE

Upholstery : Cleaning with soapy water only (no solvents), can be disinfected with products such as Phagosept Spray, Salvanios pH10 at 0.75% and Anios Surfaces Hautes, without abrasive material

Base : cleaning with soapy water only (no solvents) without abrasive materials

WARRANTY

- Mechanical/ electrical part: 5 years
- Upholstery: 1 year

CONTACT

Parc d'entreprises Visionis
278 rue de l'industrie 01090
GUÉREINS- FRANCE
COMMERCIAL SERVICE
Tél : +33(0)4 74 66 69 49

Fax : +33(0)4 74 66 69 50
export@carinamedical.com